## HiViral<sup>TM</sup> Transport Kit Art.# AL167

## **Content & Packaging**

50 or 100 round-bottom vacutainer tubes (height incl. cap: 8 cm) x 1,3 ml/ box

400 self standing tubes (height incl. cap: 10 cm) x 1,3 ml/ box

2 swabs (buccal and nasal) with predetermined breaking points

*Different tubes and swabs are available on request:* 

Tel. +49 (6251) 989 24 0 webshop@neofroxx.com

## **Storage & Shelf Life**

Store at 15-30°C before use

Transport within 24 hours at 2-8 °C

For long term storage freeze at -20 °C

(The shelf life at room temperature is 12 months)

## Composition

The HiViral<sup>™</sup> Transport medium from HiMedia is a specially formulated medium for collection, transport and long term freeze storage of viruses - designed to maintain the optimum viability and virulence of the viral sample.

HiViralTM is based on Hanks 'Balanced Salt Solution and also contains:

🖌 a protective protein,

antibiotics to prevent contamination



Phenol red is used as a pH indicator

 and a cryoprotectant which helps in preserving the viruses, if specimens are frozen for prolonged storage

